

INO Multi-Gel

Light Curing Resin for Implant Jig and Pattern



INO Multi-Gel is a gel type substance containing resin material to make implant jig, implant verification bar, pattern resin and denture resin.


INO Multi-Gel is a injection method, instead of mixing resin powder and liquid.

Usage

- Connection of support for implant prosthetics, denture working
- Inlay, Onlay, Crown Bridge Joint Connection and Making Coping and Caps
- Fabrication of partial clasp and Structures without Replication Models
- The Black out Technique

Advantages

- Easy to build up & Simple to use
- Reduce material savings and work time
- It does not need separation and polishing
- It does not have resin odor
- Easy work with moderate flow and visibility
- Less than 1% Shrinkage rate




Configuration

- MG 14 RED : 1 syringe (14ml) x 1, tips
- MG 14 RED 5 : 1 syringes (14ml) x 5, tips




Technical Data

• Depth of Curing




Light curing 30sec




Chamber curing 30sec


• Compressive yield point




• Shrinkage rate



• Flexural strength



• Elongation



INOMulti-Gel

Light Curing Resin for Implant Jig and Pattern

INOD (주)이노디
INNO DENTAL CO., LTD.

TEL : +82-2-757-1112 **FAX :** +82-2-757-9990 **E-MAIL :** info@inod.co.kr **WEB :** www.inod.co.kr
ADDRESS : #812, Byucksan Digital Valley VI, 219, Gasandigital 1-ro, Geumcheon-gu, Seoul, 08501, Republic of Korea