



INOD

U.P. FIBER SPLINT
Reinforcement UHMWPE Splint

INOD

812, Byucksan Digital Valley VI, 219,
Gasandigital 1-ro, Geumcheon-gu,
Seoul, 08501, Korea
INNO DENTAL CO., LTD.

T : +82 2 757 1112
Fax : +82 2 757 9990
E-mail : info@inod.co.kr
<http://www.inod.co.kr>

U.P. FIBER SPLINT

Reinforcement UHMWPE Splint

Clinical applications

 Orthodontic retainer	 Splinting	 Bridges	 Implant temporary
			
U.P. Fiber Splint #1 - Aesthetic - Metal free-comfortable for patients - A quick application - Easy to adapt directly to teeth and bond in place	U.P. Fiber Splint #2 - Is Quick, Clean and Simple - Is comfortable for patient (smooth and highly esthetic) - Supports teeth effectively - Has superior durability - Is easy to remove when necessary	U.P. Fiber Splint #3 - Provide strength, durability - Provide immediate convenience - Cost effective and reliable - Aesthetic - Utilize extracted teeth	U.P. Fiber Splint #4 - Aesthetic - Easy to make implant temporary prostheses - Easy to remove when necessary

Products line-up

 U.P. Fiber Splint #1 Refill 1mm x 200mm x 2pcs	 U.P. Fiber Splint #2 Refill 2mm x 200mm x 2pcs
 U.P. Fiber Splint #3 Refill 3mm x 200mm x 2pcs	 U.P. Fiber Splint #4 Refill 4mm x 200mm x 2pcs

U.P. Fiber Splint KIT

- 1mm x 200mm x 1pcs
- 2mm x 200mm x 1pcs
- 3mm x 200mm x 1pcs
- 4mm x 200mm x 1pcs
- Cutting Scissor
- Foil




U.P. Fiber Splint provides solution for periodontally compromised mobile tooth, post-trauma mobility and orthodontic treatment stabilization with non-invasive esthetic method.

Advantages

- Biocompatible splinting and orthodontic retainer
- Highly esthetic and comfortable for patient
- Very easy and fast application
- Non-invasive application requiring minimum tooth preparation
- No shelf life, no special conservation required
- Various sizes of width (1,2,3 & 4mm)
- Strong flexural strength of 380Mpa
- Excellent wettability of bonding agent
- Great bonding strength with flowable composite resin
- Affordable



The comparative graph of flexural strength for splinting products



U.P. FIBER SPLINT PLUS

Reinforcement UHMWPE Splint



This product has unfilled bonding treatment in the splint.
Remove the product and attach it directly with light curing

Products line-up

No.	Model	Length (mm)	Width (mm)	Thickness (mm)	PCS
1	U.P. FIBER SPLINT PLUS 1mm	80	1	0.2	5
2	U.P. FIBER SPLINT PLUS 2mm	80	2	0.2	5
3	U.P. FIBER SPLINT PLUS 3mm	80	3	0.2	5
4	U.P. FIBER SPLINT PLUS 4mm	80	4	0.2	5