

The masterpiece that rules all the others

Q.P. FIBER POST **PLUS**


3 reasons that you should choose Q.P. Fiber Post PLUS:

1. Excellent radiopacity and light penetration
2. High flexural strength
3. High quality with resonable price


Advantages

- Good bonding strength caused by SILANE finishing on post surface
- Higher transparency provides superior aesthetic and light transmission, and is advantageous for dual cure cementation
- Being anticorrosive and having passed all biosafety tests, it ensures bioaffinity and safety
- Passed all the tests including cytotoxicity test, ames test, acute system toxicity test and oral mucosal irritation test
- Can be atraumatically removed in minutes
- Tapered apical shape conserves dentin
- 10 different sizes and dedicated drills are available
- Elastic bending strength similar to that of dentin helps to protect teeth without causing root canal fracture
- Enhanced flexural strength ($>1,800\text{Mpa}$)
- Optimal Radiopacity ($>330\%$ AI equivalent) (ISO4049)


Technical data




Elastic modulus similar to that of dentin helps to protect teeth without causing root canal fracture



Heightened flexural strength ($>1,800\text{Mpa}$)



Optimal radiopacity ($>330\%$ AI equivalent) (ISO4049)



Excellent radiopacity



Perfect light penetration



Safe drill of reamer shape

Products line-up



Description KIT Options

Kit	Diam	Q.P Fiber Post PLUS	Drills
QPA1209	ø1.14	10pcs	1pcs
	ø1.25	10pcs	1pcs
	ø1.55	10pcs	1pcs

Kit	Diam	Q.P Fiber Post PLUS	Drills
QPA1309	ø0.9	5pcs	1pcs
	ø1.0	5pcs	1pcs
	ø1.14	5pcs	1pcs
	ø1.25	5pcs	1pcs


Kit	Diam	Q.P Fiber Post PLUS	Drills
QPA1409	ø1.0	5pcs	1pcs
	ø1.2	5pcs	1pcs
	ø1.4	5pcs	1pcs
	ø1.55	5pcs	1pcs

Description of Refills

Model	QPA9	QPA10	QPA14	QPA20	QPA25	QPA40	QPA50	QPA55	QPA80	QPA200
Size	ø0.9	ø1.0	ø1.14	ø1.20	ø1.25	ø1.40	ø1.50	ø1.55	ø1.8	ø2.0
Refill	10pcs	10pcs	10pcs	10pcs	10pcs	10pcs	10pcs	10pcs	10pcs	10pcs
Combination with Drill										
Removal of the Root Filling	Universal Drill PDO									
Preparation of the Root Canal	PD9	PD10	PD1	PD20	PD2	PD40	PD50	PD3	PD80	PD200
Color Code	Brown	Yellow	Blue	Green	Red	Purple	Black	Orange	White	Gray



Diameters in millimeters



INOD

INNO DENTAL CO., LTD.

812, Byucksan Digital Valley VI, 219, Gasan-dong 1-ro,
Geumcheon-gu, Seoul, 08501, Korea

<http://www.inod.co.kr>

T : +82 2 757 1112 / Fax : +82 2 757 9990 / E-mail : info@inod.co.kr